






PICO sNPWT

Product list by SKU

Dressings	Dressing Sizes	PICO 7 System		PICO 14 System	Multipack with	PICO 7Y device
		+ 1 dressing	+ 2 dressings	+ 2 dressings	5 dressings	+ 2 dressings
		Code	Code	Code	Code	Code
	20cm x 25cm	-	-	-	-	66022031
	10cm x 20cm	66022012	66022002	66022042	66022022	-
	10cm x 30cm	66022013	66022003	66022043	66022023	-
	10cm x 40cm	66022014	66022004	66022044	66022024	-
	15cm x 15cm	66022015	66022005	66022045	66022025	-
	15cm x 20cm	66022016	66022006	66022046	66022026	-
	15cm x 30cm	66022017	66022007	66022047	66022027	-
	20cm x 20cm	66022018	66022008	66022048	66022028	-
	25cm x 25cm	66022019	66022009	66022049	66022029	-

Advanced Wound Management
Smith & Nephew
Fort Worth, TX 76109
USA
T +1 (800) 876 1261
F +1 (727) 392 6914

www.smith-nephew.com
www.possiblewithpico.com

◇ Trademark of Smith+Nephew
©2022 Smith+Nephew
PCPE12-34134-0322

Important Safety Information

The PICO pumps contain a MAGNET. Keep the PICO pumps at least 4 inches (10 cm) away from other medical devices at all times. As with all electrical medical equipment, failure to maintain appropriate distance may disrupt the operation of nearby medical devices. For full product and safety information, please see the Instructions for Use.